

dexcomG6

MAKE KNOWLEDGE YOUR SUPERPOWER WITH DEXCOM G6



Dexcom G6 receiver (not shown) is part of the Dexcom G6 CGM System. Smart devices sold separately.
To view a list of compatible smart devices, visit dexcom.com/compatibility.

“I’M A BIG KID AT HEART AND DEXCOM HELPS ME LIVE MY LIFE!”

Brandon D.
Dexcom Warrior with diabetes

EMPOWERED DIABETES MANAGEMENT STARTS HERE.

Dexcom G6 provides real-time glucose readings as often as every five minutes. It connects with your compatible smart device* to keep you aware of where your glucose is now and where it’s heading, whether you’re at home or on the go.



Zero Fingersticks

Factory calibrated to allow for zero fingersticks†

†If your glucose alerts and readings from Dexcom G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



Customizable Alerts

Set your range and get notified when you’re high or low



Slim, Water-Resistant Sensor†

A slim, water-resistant sensor is discreet and easy to wear for up to 10 days



Indicated for Ages 2+

The only CGM FDA-permitted for children ages 2 years and up



A1C

Keeping You in Control

Dexcom CGM use can help you reduce A1C and improve glycemic control^{1,2}



Simplified Setup

Easy-to-understand product instructions and painless sensor insertion³

TAKE A CLOSER LOOK AT THE DEXCOM G6 CONTINUOUS GLUCOSE MONITORING SYSTEM



Smart devices sold separately.

“NO FINGERSTICKS.† BEST. NEWS. EVER.”

Dexcom G6 eliminates fingersticks for calibration and diabetes treatment decisions.†

†If your glucose alerts and readings from Dexcom G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



Anita V.
Dexcom Warrior with diabetes

*To view a list of compatible smart devices, visit www.dexcom.com/compatibility. †Dexcom G6 sensors and transmitters are water-resistant and may be submerged under eight feet of water for up to 24 hours without failure when properly installed. 1 Beck RW et al. JAMA 2017;317:371-8. 2 Welsh JB, et al. Diabetes Technol Ther 2019;21:128-32. 3 Shah VN, et al. Diabetes Technology and Therapeutics. 2018;20(6): 428-433.

PROVEN ACCURACY¹ PUTS YOU BACK IN CONTROL

Dexcom G6 gives you the flexibility and convenience to monitor your glucose data when, where and how you want.

HOW DOES THE DEXCOM G6 CGM SYSTEM WORK?



1 APPLICATOR

- Simple auto-applicator inserts the sensor just beneath the skin

2 SENSOR

- Measures glucose levels for up to 10 days
- Acetaminophen blocking: provides accurate glucose readings, even when taking acetaminophen^{* 2}

3 TRANSMITTER

- Affixed to the top of sensor; sends data wirelessly to a display device

4 DISPLAY DEVICE

Dexcom receiver

- Touchscreen display
- Displays a glucose reading as often as every five minutes — up to 288 a day
- Trend arrows and vivid colors show glucose speed and direction your glucose is heading

Compatible Smart Device[†] (Not included with Dexcom G6 CGM System)

- Check your glucose values from work, home, or school using the Dexcom G6 app on a compatible smart device

STAY MOBILE WITH THE DEXCOM G6 APP

Manage your diabetes from the palm of your hand. The Dexcom G6 app allows you to quickly access your glucose information, while the customizable alerts can notify you of highs and lows.



The Dexcom G6 app is compatible with both iOS and Android devices.[†]



^{*}Dexcom G6 readings can be used to make diabetes treatment decisions when taking up to a maximum acetaminophen dose of 1,000mg every 6 hours. Taking a higher dose may affect the Dexcom G6 readings. [†]To view a list of compatible smart devices, visit dexcom.com/compatibility. [§]Separate Follow App required. 1 Shah VN, et al. Diabetes Technology and Therapeutics. 2018;20(6): 428-433. 2 Calhoun P, et al. J Diabetes Sci Technol 2018;12:393-6.

THE DEXCOM SHARE FEATURE LETS YOU TAKE THE LEAD, SO OTHERS CAN FOLLOW

The built-in Share feature of the Dexcom G6 app allows up to 10 people to follow your glucose levels, giving you the circle of support you need.* By downloading the Dexcom Follow¹ app, followers can view your glucose data directly from their smart device, whether they live down the street or across the country.



Zola & Keary C.
Dexcom Warriors
with diabetes



Smart devices sold separately.†

*Do not make treatment decisions based upon Dexcom Share and Dexcom Follow readings. Always confirm with your compatible smart device or receiver. †Internet connectivity required to access Dexcom Follow. ‡To view a list of compatible smart devices, visit dexcom.com/compatibility.

TRACK YOUR PROGRESS WITH DEXCOM CLARITY

Dexcom CLARITY puts the spotlight on your glucose insights, helping to inform you and your diabetes care team. It allows you to share data with your clinic and monitor glucose patterns, trends, and statistics in between visits. See what's possible with Dexcom CLARITY!

WHEN AT HOME



Log in at clarity.dexcom.com


Use your current Dexcom login or, if needed, create an account.


WHILE ON THE GO



Get weekly notifications with the Dexcom CLARITY app

Notifications available to Dexcom mobile users.





“MY A1C HAS NEVER BEEN LOWER...THAT’S THE BEAUTY OF ALWAYS KNOWING WHERE YOU ARE.”

-Vickie G.
Dexcom Warrior
with diabetes

Individual results may vary.



“YOU HAVE GIVEN MY HUSBAND AND ME OUR LIVES BACK. WE’VE RETURNED TO WORK, ENROLLED THEM IN DAYCARE, AND NOW CAN SLEEP THROUGH THE NIGHT!”

-Rebeka J.
mom of daughters with
diabetes, Dexcom Warriors



“MY DEXCOM G6 KEEPS ME IN THE GAME!”

-Nathan P.
Dexcom Warrior
with diabetes

“DEXCOM SHARE IS LIFE-CHANGING AS SHE GOES OFF TO SCHOOL!”

-Beth R.
mom of daughter with
diabetes, Dexcom Warriors



“THE DEXCOM G6 HAS PROVEN TO BE AN INCREDIBLY USEFUL TOOL IN MY DIABETES CARE AND MANAGEMENT.”

-Molly J.
Dexcom Warrior
with diabetes



“NOT HAVING TO PRICK MY FINGERS FEELS AMAZING!”

-Lala J.
Dexcom Warrior
with diabetes

If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



“MY WHOLE LIFE IS ON MY PHONE... SO THE DEXCOM G6 FITS RIGHT IN.”

-Harrel C.
Dexcom Warrior
with diabetes



“I REALLY DON’T UNDERSTAND HOW YOU CAN CONTROL DIABETES WITHOUT A DEXCOM SYSTEM.”

-Gary K.
Dexcom Warrior
with diabetes



Dexcom
WARRIORS

DEXCOM SUPPORT

From paperwork to insurance processing, to product support and training, our Customer Service and Patient Support teams are here to assist you every step of the way.



DEXCOM CARE

Dexcom CGM training, software downloads, and tutorials

1-877-339-2664

Monday - Friday 6 AM - 7 PM PST (subject to change)

Saturday 6 AM - 2:30 PM PST (subject to change)

patientcare@dexcom.com



CUSTOMER SALES SUPPORT

Support with Dexcom orders and general customer questions

1-888-738-3646

Monday - Friday 6 AM - 5 PM PST

Saturday 8 AM - 2 PM PST

customerservice@dexcom.com



TECHNICAL SUPPORT

Product troubleshooting or replacement inquiries

1-844-607-8398

Available 24 hours a day; 7 days a week

techsupport@dexcom.com

Visit dexcom.com/G6 for more information.

NEXT STEPS

Ready to start the process? Follow these simple steps.

1 FILL OUT THE FORM

Please fill out the form included in your brochure and give it to your doctor's office staff to start the process to get your Dexcom CGM.

2 LEARN MORE ABOUT DEXCOM G6

Visit dexcom.com/G6 to learn more about the benefits of Dexcom CGM.

3 TALK TO DEXCOM

A representative from Dexcom will call you at the number provided on your form. It may come up on your phone as a local number. Please listen to the message and stay on the line to speak with a Dexcom representative.

4 GET A FREE BENEFITS CHECK

Dexcom will conduct a free benefits check on your behalf and relay your out-of-pocket costs before we ship. Once you agree to your costs, we will arrange for shipment.

5 START YOUR DEXCOM

Once your boxes arrive, open your Dexcom G6 CGM System and follow the instructions. We also have video tutorials online at dexcom.com/welcome

If you have any questions, please call Dexcom CARE and speak with a Diabetes Coach at **1-877-339-2664**.

Dexcom

BRIEF SAFETY STATEMENT

Failure to use the Dexcom G6 Continuous Glucose Monitoring System (Dexcom G6) and its components according to the instructions for use provided with your device and available at <https://www.dexcom.com/safety-information> and to properly consider all indications, contraindications, warnings, precautions, and cautions in those instructions for use may result in you missing a severe hypoglycemia (low blood glucose) or hyperglycemia (high blood glucose) occurrence and/or making a treatment decision that may result in injury. If your glucose alerts and readings from Dexcom G6 do not match symptoms or expectations or you're taking over the recommended maximum dosage amount of 1000mg of acetaminophen every 6 hours, use a blood glucose meter to make diabetes treatment decisions. Seek medical advice and attention when appropriate, including for any medical emergency.